



### FlashMax P4 Light Curing Pen

- Highest output intensity of LED curing lights
  - Light tip can adjust  $\pm 15^\circ$  for better access to distal regions of the mouth; 8 mm light tip is ideal for use with fixed appliances.
  - Economical unit provides less heat without sacrificing power and efficiency
  - Treat the upper and lower arches with maximum efficiency using the bilateral control buttons. Identical control buttons on either side allow you to rotate the light to achieve the best angle for maximum bond strength
  - A battery-operated, cordless hand unit that provides maximum office flexibility and convenience
  - Three curing modes for enhanced bonding versatility and operating efficiency
- Peak Output Intensity: 6,000 mW/cm<sup>2</sup>
  - Curing Modes: 1s/3s
  - Wavelength: 440 nm – 480 nm
  - Weight : 120g
  - Battery: LiFePO<sub>4</sub>
  - Recharge Time: 2.5 hours
  - Dimensions: 225mm x 25 mm
  - Manufacturer Warranty: 2 years

Description	Package Contains	Order Number
FlashMax P4 Light Curing Pen Ortho PRO	1 each	J02130

### FlashMax Accessories



Description	Package Contains	Order Number
4 mm tips	50 each	J02121
8 mm tips	50 each	J02122
Covers	10 each	J02123
Shield	1 each	J02124
Universal disposable light shield	50 each	J02125
Disposable light sleeves	250 each	J02126

### Mini Led Ortho 2

Cure in a simple gesture, in few seconds from few brackets to a whole quadrant

- Large photo-initiator compatibility**  
 Relevant wavelength spectrum: Camphoroquinone (CQ) 470nm but also PPD and PAB 430 nm  
 Large power range from 1250mW/cm<sup>2</sup> to 3000mW/cm<sup>2</sup>  
 Customizable modes to polymerize diverse materials in various clinical situations
- Ortho friendly**  
 Customizable exposure times of 4, 8, 12 and 32 sec, specifically studied to move from one bracket to another or to cure a full arch
- Customizable signal**  
 An alert (beep, micro-flash, mixed or no alert) prompts the user to move to the next curing position every 4 sec
- Well-balanced combination**  
 Cordless light handpiece easy to grip  
 Robust body in aluminium  
 Li-Ion battery with no memory effect
- Efficacy and safety**  
 Cure any type of material used in orthodontics efficiently and quickly. The absence of UV reduces heat production for maximum comfort and safety
- Long lasting battery**  
 MiniLED Ortho 2 can cure 15 patients at maximum power, no need to recharge it during the day



Description	Package Contains	Order Number
Mini Led Ortho 2 Light Pen	1	ECM8570
Opalescent Tip Ø 5,5 mm, 45°	1	ECM8565

### Light Pen Glasses

Glasses for light pen or plasma

Package Contains	Order Number
1	ECM8521

