



### Glass Ionomer Cement for Orthodontic Bonding GC Fuji

#### GC Fuji Ortho

- Excellent bond strength even in a humid environment: easy application
- Easy removal of orthodontic appliances: no enamel fractures
- Self-curing
- Long term fluoride release
- No discoloration

Description	Package Contains	Order Number
Powder 40 gr/liquid 2 x 6,8 ml	1	ECM8084
Liquid 23,8 ml	1	ECM8087
Powder 40 gr (refill)	1	ECM8093
Liquid 6,8 ml (refill)	1	ECM8094



#### GC Fuji Ortho LC

- Light curing within 40 sec
- Excellent adhesion even in a humid environment
- Easy removal of orthodontic appliances
- Long term fluoride release
- No discoloration

Description	Package Contains	Order Number
Capsules	50	ECM8085
Powder 15 gr / Liquid 6,8 ml	1	ECM8086
Powder 40 gr (refill)	1	ECM8095
Liquid 6,8 ml (refill)	1	ECM8096



#### GC Fuji OrthoBand

- For bonding of bands and appliances
- Paste pak cartridge delivery
- Optimal handling
- High fluoride release
- Easy removal during application and after debonding

Description	Package Contains	Order Number
Paste pak (2 cartridges 13,3 fr + mixing pads)	1	ECM8092

## Glass Ionomer Cement for Orthodontic Bonding GC Fuji

### GC Fuji I

The biocompatible, resin-free luting solution with 25 years of great clinical success. Fuji I, GC's first luting GI cement, is widely used across the world.

- Chemical adhesion to tooth structure, ensuring an excellent marginal seal
- A moisture-tolerant solution when a dry field cannot be assured
- Quick and easy without the need for etching and bonding of the tooth
- Rubbery consistency during setting for an easy removal of excesses
- Resin-free cement for patients allergic to methacrylates
- Virtually no post-operative sensitivity: a pulp-friendly solution



Description	Package Contains	Order Number
(Powder 15 g, liquid 7 ml, conditioner 6,5 ml & accessories)	1	ECM8088
Package of 50 capsules	1	ECM8090

### GC Fuji PLUS

Fuji PLUS has been designed to offer you simple and error-free cementation procedures and a strong and durable bond you can rely on. The outstanding performance and the availability in various shades make GC Fuji PLUS a multi-purpose cement.

- A moisture-tolerant solution when a dry field cannot be assured
- Quick and easy without a need for complicated bonding procedures
- Rubbery consistency for easy removal of excess
- Pulp-friendly with virtually no post-operative sensitivity
- Proven and predictable performance for most routine cases
- Also available in extended working time (EWT) version for long-span bridges



Description	Package Contains	Order Number
Package of 50 capsules	1	ECM8081
Full kit (powder 15 g, liquid 7 ml, conditioner 6,5 ml & accessories)	1	ECM8082
Capsule applicator	1	ECM8089
Capsule mixer	1	ECM8097

## Glass Ionomer Luting Cement : Aqua Cem

For bands fixation, powder 30g.

Package Contains	Order Number
1	ECM8013

